

H1N1 Probe List

Ready-made RT-PCR probe, designed fully based on WHO-confirmed sequence and used in influenza A (H1N1) detection area, now is available for your work at any time!

InfA-F	GACCRATCCTGTCACCTCTGAC
InfA-R	AGGGCATTYTGGACAAAKCGTCTA
InfA-probe	FAM-TGCAGTCCTCGCTCACTGGGCACG-BHQ1
SWInfA-F	GCACGGTCAGCACTTATYCTRAG
SWInfA-R	GTGRGCTGGGTTTTTCATTTGGTC
SWInfA-probe	FAM-CYACTGCAAGCCCAT(BHQ1)ACACACAAGCAGGCA-PO4
SWH1-F	GTGCTATAAACACCAGCCTYCCA
SWH1-R	CGGGATATTCCTTAATCCTGTRGC
SWH1-probe	FAM-CAGAATATACAT(BHQ1)CCRGTCACAATTGGARAA-PO4
RnaseP-F	AGATTTGGACCTGCGAGCG
RnaseP-R	GAGCGGCTGTCTCCACAAGT
RnaseP-probe	FAM-TTCTGACCTGAAGGCTCTGCGCG-BHQ1



Shanghai ShineGene Molecular Bio-tech Co.,Ltd.

Add: Floor 2, Building A, 328#, Wuhe Road, Shanghai 201109,

Tel: +86-21-54460832

Fax: +86-21-54460831

E-mail: shinegene@vip.163.com

Website: www.synthesisgene.com