

Certificate of dATP Analysis

Product :	2'-Deoxyadenosine -5'-Triphosphate, Sodium Salt (dATP) Sequenase and PCR Qualified
Storage Recommendations:	Always avoid freeze-thaw cycles or exposure to frequent temperature changes. These fluctuations can greatly alter product stability.
Long term (Infrequent use: 1-2 times per month):	-70°C
Daily/Weekly use:	-20°C
Molecular Formula:	$C_{10}H_{12}N_5O_{12}P_3Na_4$
Molecular Weight:	579.2
Form:	Aqueous solution, 100mM
Purity (HPLC):	>99.8%
λmax:	259nm ± 2nm
Absorbance Rations (pH 7.0):	A250/A260 0.80±0.03
	A280/A260 0.15±0.03
PH:	7.5 ± 0.1
Quality control:	The concentration is verified by optical density spectrophotometry. This preparation is free of DNase and RNase contamination as
Shine Gene	determined by incubation of the 4 dNTPs with radioactive substrates.
	3. The purity of this preparation is determined by HPLC. The minimum passing specification for the triphosphate content is 99.8%.
	4. This lot of dNTPs has been performance tested with Sequenase, <i>Taq</i> DNA polymerase, <i>Pfu</i> DNA polymerase and <i>T7</i> DNA polymerase.